



Cleaning Procedure for Gas Sampling Bags

Adtech gas sample bags can be used on multiple occasions, however, they should be cleaned between each sampling cycle, to prevent contamination.

Once the testing cycle on the fluoroplastic bag has been completed, the bags need to be flushed through (purged) with clean air.

This is a simple process:

- ◆ Remove the air that was being sampled from the bag.
- ◆ Fill the bag with clean air.
- ◆ Hold the clean air within the bag for 2-3 minutes.
- ◆ Remove the air from the bag.
- ◆ Carry out a leak test on the bag.

This process can be carried out as many times as required.

Most analysing equipment will have the ability to purge the bags and carry out a leak test.

This information is a guideline only and you may wish to modify the cleaning process to suit each application.

For more information about our Gas Sample Bags, please visit [adtech.co.uk](https://www.adtech.co.uk).

For technical assistance, contact our customer service team:

+44 (0)1285 762000 - sales@adtech.co.uk